

Introduction

In monoclonal antibody (mAb) manufacturing, accurate antibody titre determination and monitoring are crucial for clonal selection, optimisation of expression rates, process control and scale-up. For accurate antibody titre quantification, Protein-A Affinity HPLC offers excellent reproducibility, robustness and precision over ELISA-based assays.

The SOLAD[™] Protein-A is a novel affinity HPLC column, based on the next-gen SOLAD[™] platform technology, for quick separation and robust quantification of various antibodies. It provides rapid, accurate titre analysis with linearity across a wide concentration range, without any carryover.

SEM of SOLAD Particles (2500x magnified)



Highlights

- Wide dynamic range
- Fast Analysis
- High sensitivity & specificity for mAbs
- ✤ Long lifetime

SEM of SOLAD Particles (10000x magnified)



Column Technology

SOLAD[™] is a proprietary product range of highly monodisperse, spherical, durable and reproducible non-porous silica particles manufactured by Glantreo. For protein and antibody analysis, SOLAD[™] offers excellent mechanical and chemical stability, with minimal fouling and non-specific adsorption, compared to porous silica and polymeric particles.

The SOLAD[™] Protein-A column combines 5µm and 10µm sized SOLAD particle with recombinant Protein A ligand through proprietary cutting-edge bonding technology. The result being a hydrophilic resin surface with a high pH stability and a rugged and reproducible column. The column offers superior efficiency and resolution, without analyte carryover.

The SOLAD[™] Protein-A column can accurately measure mAb titres at various phases of development due to its high loading capacity and excellent repeatability. This column is also suitable for small-scale purification of mAbs for initial characterisation due to its reduced Protein A leaching. The column is available in PEEK and Bioinert Stainless Steel column hardware in various dimensions.



Glantreo Limited ERI Building, Lee Road, Cork, Ireland, T23 XE10

Contact Details Tel: +353 21 490 1965 Email: info@glantreo.com www.glantreo.com





DISCLAIMER: The information set forth in this Product Data Sheet represents typical properties of the product described; the information and the typical values are not specifications. Glantreo Limited makes no representation or warranty concerning the products, expressed or implied, by this Product Data Sheet.

Technical Specifications and Applications of SOLAD[™] Protein-A Column



Rapid separation of IgG from impurities







Figure 3: Effect of Flow Rate on IgG binding

Column – SOLAD™ Protein-A Column, 4 x 50 mm PEEK Binding buffer (A) - 50 mM Sodium Phosphate pH 7.4 Elution buffer (B) - 50 mM Sodium Phosphate pH 2.5 Gradient – 0%B: 0-0.5 min, 100%B: 0.5-1.1 min, 0%B: 1.1-4 min

Flow rate - 2 mL/min, Detection: UV @ 280 nm, Temperature - 25 C

Sample – 20 µL cell culture supernatant containing mAb





Figure 1 depicts a 20 µL injection of mAb feedstock onto the SOLAD[™] Protein A column. The proprietary column technology yields sharp, concentrated peaks and fast, efficient elution.

As can be seen in Figure 2, SOLADTM Protein-A delivers accurate titre measurement with linearity, with a dynamic loading capacity in excess of $300 \ \mu g$.

The SOLADTM Protein-A column can be used from flow rates upto 2.5 mL/min, as seen in Figure 3. This allows fast separation of antibodies without loss of efficiency.

Features of SOLAD™ Protein-A Column	
Support Matrix	SOLAD [™] non-porous hydrophilic resin
Ligand	Recombinant Protein-A
Particle sizes	5 μm 10 μm
Monodispersity Index (d90/d10)	1.1
Flow rate	≤ 2.5 mL/min
Pressure Limit	2000 psi
Temperature	≤35 °C
pH range	2-8
IgG sample load	≤300 ug (for 4 x 50 mm column)
Column Hardware Material	PEEK, Bio-inert Stainless Steel
Column Hardware Dimension Options	PEEK:
	◆ 2.1 x 30 mm ◆ 2.1 x 50 mm ◆ 4.0 x 50 mm
	Bio-inert Stainless Steel: • 2.1 x 30 mm • 4.6 x 50 mm

Glantreo Limited

ERI Building, Lee Road, Cork, Ireland, T23 XE10

Contact Details Tel: +353 21 490 1965 Email: info@glantreo.com www.glantreo.com





DISCLAIMER: The information set forth in this Product Data Sheet represents typical properties of the product described; the information and the typical values are not specifications. Glantreo Limited makes no representation or warranty concerning the products, expressed or implied, by this Product Data Sheet.